


EPR-spectroscopy in photochemistry



The EPR technique is one of the most sensitive physical methods of detection, identification and study of free radicals, ion-radicals and triplet-state molecules generated during photochemical processes such as:


- 1) photodecomposition (homolytic photolysis into radicals, heterolytic photolysis into ions — photoionization);
- 2) electron phototransfer (photooxidation, photoreduction, photodissociative electron addition);
- 3) photosensibilization (photosensitized oxidation, photosensitized reduction, photosensitized decay) [1].



EPR spectroscopy was used to illustrate that degradation of organic Benzophenone-3 (BP-3) UV filter using $\text{Sb}_2\text{O}_3/\text{TiO}_2$ and PbO/TiO_2 photocatalysts under ultraviolet irradiation proceeds by a radical mechanism. The presence of a peak of the TEMP-1O2 adduct on the spectra indicates the participation of singlet oxygen in the degradation of the UV filter in sunscreen, which was detected in surface and groundwater [2].

Polymer materials for photovoltaics are subject to photodegradation, which occurs under the influence of UV radiation and leads to a decrease in conductive properties. Using the PS 100X EPR spectrometer and irradiation of the test sample with visible light directly in the resonator of the device in the temperature range 90-340K, the nature and environment of paramagnetic centers arising in a sample of a composite fullerene-containing polymer were studied [3].

EPR-based research [4] on CMS 8400 spectrometer showed, that the hydroxyl radical ($\cdot\text{OH}$), generated during the photo-Fenton reaction mediated by hybrid materials catalyst $\alpha\text{-Fe}_2\text{O}_3\text{-Graphene}$, is a major reaction product responsible for the rhodamine degradation – a model wastewaters pollutant.







Figure 1 – EPR-spectra of UV-generated DMPO-OH (a) и DEPMPO-OH (b) radicals.

Experiment parameters: center field 336.73mT; sweep width 12mT; modulation frequency 109.375kHz; modulation amplitude 200uT; attenuation 25dB; number of points 1000; sweep time 60s.

REFERENCES

1. Kholmogorov V. E. THE METHOD OF ELECTRON PARAMAGNETIC RESONANCE IN PHOTOCHEMISTRY. // ADVANCES IN CHEMISTRY – 1968– T. XXXVII– vol.8 –С. 1492-1520;
2. Z. Wang, A. Deb, V. Srivastava, S. Iftekhar, I. Ambat, M. Sillanpää. - Investigation of textural properties and photocatalytic activity of PbO/TiO₂ and Sb₂O₃/TiO₂ towards the photocatalytic degradation Benzophenone-3 UV filter. - Separation and Purification Technology, 228 (2019) 115763-115775
3. V. I. Krinichnyi and E. I. Yudanova Light-induced EPR spectroscopy of charge carriers photoinduced in polymer/ fullerene bulk heterojunctions JOURNAL OF RENEWABLE AND SUSTAINABLE ENERGY 1, 043110_2009_18p.
4. S. Frindy, M. Sillanpää Synthesis and application of novel α -Fe₂O₃/graphene for visible-light enhanced photocatalytic degradation of RhB // Materials & Design–2020–Vol. 188, 108461

To demonstrate the application of EPR technique in photochemistry research the hydroxyl radicals UV-generation process has been used. The radicals have been detected on bench-top EPR spectrometer Spinscan X. Spin-trapping is the method for short-living radicals (such as the OH-radicals) detection. DMPO (3,4-dihydro-2,3-dimethyl-2H-pyrrole 1-oxide) и DEPMPO (P-(3,4-dihydro-2-methyl-1-oxido-2H-pyrrol-2-yl)-phosphonic acid, diethyl ester) solutions have been used as spin traps.

Hydroxyl radicals are generated in the reaction mixture consisting of 10mM H₂O₂ with 100mM DMPO or 50mM DEPMPO, under the UV irradiation with $\lambda=395$ nm. The DEUTERIUM TUNGSTEN HPOWER/LS-DWHP lamp was used as a UV-source, the light was directed onto the sample by using the optical fiber cable through the EPR-resonator window.

Numerous examples of EPR-based experimental research of photochemical processes demonstrate the broad scope of practical applications of the EPR technique for studying the irradiated products. The EPR spectroscopy can be used for deep scientific explorations in photochemistry and spin-photonics as well as for the development of spin-based technologies, particularly: photolysis control, photonics and spintronics, - the frontier research areas allowing for new methods of photoinduced processes research.



Our solution –
EPR spectrometer
Spinscan X